



# CARDIOS



# about Cardios

*Founded in 1977, Cardios offers products and services to diagnose and prevent heart and blood vessel related diseases, like hypertension, heart attacks, and heart arrhythmia. Our R&D team develops state-of-art digital Holter and blood pressure monitors and also systems for analysis and telemedicine connection, that are developed in our own certified production line. After sales support is provided to customers and includes scientific and technical assistance and service repair center.*

*Cardios' headquarters and operations are based in Brazil. As part of Halma's Medical Sector, Cardios' R&D team brings innovation and new technologies, opening new opportunities of international expansion and business models.*



# recorder

*Brings innovation to Holter by using the latest ECG technology in a modern, compact form*

# Cardio Light+

**24-72h Digital Holter Recorder**



A signal acquisition rate of 800 samples per second, pace maker detection and a real time internal clock. Continuous recording for up to 72h without battery replacement.



Data can be transmitted over the internet for remote analysis, case discussions or database storage through CardioNet.



Its reduced size and weight, along with its award winning design, makes it easily portable and comfortable to patients.



## Technical Specifications

Size	82 x 60 x 14mm
Weight	47 g
Storage	Internal memory (USB-C)
Power supply	1 rechargeable battery (AAA)
Signal acquisition	800sps - 12 bits

**€ 1434**



Av. Paulista, 509 1º andar  
01311-910 São Paulo SP Brazil

+55 11 3883-3000  
+55 11 3883-3030  
[sales@cardios.com.br](mailto:sales@cardios.com.br)  
[www.cardios.com.br](http://www.cardios.com.br)



**Obelis s.a.**  
Bd. Général Wahis 53  
1030 Brussels, Belgium

+ 32 2 732 59 54  
+ 32 2 732 60 03  
[mail@obelis.net](mailto:mail@obelis.net)



A **Halma** company